



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

L-HISTIDINE MONOHYDROCHLORIDE MONOHYDRATE TESTING METHODS

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

Page **1** of **12**

Controlled Copy Number: 2, Controlled Copy Location: Website, Printed By: CATHERINE.SUTTON, on 13 Jan 2026

TABLE OF CONTENTS

1.	PURPOSE:.....	3
2.	SCOPE:.....	3
3.	RESPONSIBILITIES:	3
4.	EQUIPMENT:	3
5.	REAGENTS:	3
6.	REFERENCES:	4
7.	ANALYTICAL PROCEDURES:.....	5
	TABLE 1: CHROMATOGRAPHIC SYSTEM.....	8
	TABLE 2: CHROMATOGRAPHIC SYSTEM.....	9
8.	COMPENDIAL DIFFERENTIATIONS:	12
	TABLE 3: COMPENDIAL ANALYSES NOT HARMONIZED	12
	TABLE 4: ONE METHOD NAME.....	12
	TABLE 5: IN-HOUSE VALIDATED METHODS.....	12
	TABLE 6: IN-HOUSE METHODS.....	12
	TABLE 7: OUTSIDE APPROVED LABORATORY TESTING.....	12

1. PURPOSE:

- 1.1. To provide Laboratory Personnel with procedures for testing L-Histidine Monohydrochloride Monohydrate.

2. SCOPE:

- 2.1. Applies to the testing of L-Histidine Monohydrochloride Monohydrate in the laboratory at all BioSpectra Facilities. Methods include testing for all types of L-Histidine Monohydrochloride Monohydrate sold by BioSpectra; only the specific tests required for the desired type must be tested.

3. RESPONSIBILITIES:

- 3.1. The Laboratory Manager is responsible for the control, training, maintenance, and implementation of this procedure.
- 3.2. The Laboratory Technicians are responsible for compliance with the terms of this procedure. This includes notifying the Quality Assurance / Laboratory Managers or designee if any analyses fail to meet their respective specifications.
- 3.3. Standard laboratory safety regulations apply. Before working with any chemical, read and understand the Safety Data Sheet (SDS).

4. EQUIPMENT:

- 4.1. Analytical Balance
- 4.2. Anton Paar MCP 5300 Polarimeter
- 4.3. Calibrated Pipettes
- 4.4. Hach Portable Turbidimeter
- 4.5. Metrohm Titrando 907 Auto-Titrator
- 4.6. pH Probe
- 4.7. XL200 pH/mV/Conductivity Meter, or equivalent
- 4.8. OPI-180 OD Handheld Colorimeter

5. REAGENTS:

- 5.1. **Acetic acid:** Purchased Commercially.
- 5.2. **Ammonia Water:** See Ammonia TS.
- 5.3. **Ammonia TS:** Purchased Commercially.
- 5.4. **0.1N Sodium Hydroxide:** Purchased Commercially.
- 5.5. **1M Sodium Hydroxide:** Purchased Commercially.
- 5.6. **0.1N Silver Nitrate:** Purchased Commercially.
- 5.7. **1-Propanol:** Purchased Commercially.
- 5.8. **Acetic Acid, Dilute:** 12g of glacial acetic acid diluted to 100mL with water.
- 5.9. **Barium Chloride TS:** Dissolve 12 grams of Barium Chloride Dihydrate in purified water. Filter and dilute to make a total volume of 100mL with purified water.
- 5.10. **Butanol:** Purchased Commercially.
- 5.11. **Dilute Nitric Acid:** Dilute 20 grams of Concentrated Nitric Acid to 100mL with purified water.
- 5.12. **Dipotassium Sulfate:** Purchased Commercially.
- 5.13. **EP Reference Solution BY₆:** Purchased Commercially.
- 5.14. **Hydrochloric Acid, concentrated:** Purchased Commercially.
- 5.15. **Hydrochloric Acid R1:** Dilute 70 grams of concentrated Hydrochloric Acid to 100mL with Purified Water and mix well.
- 5.16. **Hydrochloric Acid, 6N:** Pipette 51.5mL of concentrated hydrochloric acid into a 100mL volumetric flask containing a small amount of purified water. Dilute to volume with purified water.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 5.17. **Hydrochloric Acid, dilute (10%):** Dilute 23.6mL of hydrochloric acid with purified water to make 100mL.
- 5.18. **LAL Reagent Water:** Purchased Commercially.
- 5.19. **L-Histidine Monohydrochloride Monohydrate CRS:** Purchased Commercially.
- 5.20. **Ninhydrin:** Purchased Commercially.
- 5.21. **Ninhydrin Solution:** Dissolve 0.2g of ninhydrin in 5 mL of dilute acetic acid and 95 mL of butanol.
- 5.22. **Nitric acid, dilute:** Pipette 10.5mL of concentrated nitric acid into a 100ml volumetric flask containing a small amount of purified water. Dilute to volume with purified water.
- 5.23. **Potassium Hydrogen Phthalate (KHP):** Prepare an appropriate sample container at 120°C for 30 minutes. Allow to cool in a desiccator. Crush and dry a suitable amount of Potassium Hydrogen Phthalate. Dry at 120°C for 2 hours. Cool and store in a desiccator in a closed container. Stable for 3 months.
- 5.24. **Purified Water:** In-House or Purchased Commercially.
- 5.25. **Reagent Alcohol:** Purchased Commercially.
- 5.26. **Reagent Alcohol (30% V/V):** Pipette 30 mL of reagent alcohol into a 100 mL volumetric flask. Dilute to volume with purified water and mix.
- 5.27. **Sodium Hydroxide, 200g/L:** Weigh 20g of sodium hydroxide and transfer into a 100 mL volumetric flask. Dissolve and dilute to volume with purified water.
- 5.28. **Sodium Nitrite Solution:** Weigh 1g of Sodium Nitrite, transfer to a 10mL volumetric flask and dilute to volume with purified water. Prepare immediately before use.
- 5.29. **Sulfanilic Acid:** Purchased Commercially.
- 5.30. **Sulfuric Acid, concentrated:** Purchased Commercially.
- 5.31. **Sulfuric Acid, 0.005M:** Dilute 10mL of 0.05M (0.1N) Sulfuric Acid to 100mL in a volumetric flask.
- 5.32. **250g/L Barium Chloride:** Weigh 25.0 grams of barium chloride and transfer to a 100mL volumetric flask. Dissolve and dilute to volume with purified water.
- 5.33. **0.01-1.0 EU/mL Endotoxin Cartridges:** Purchased Commercially.

6. REFERENCES:

- 6.1. BSI-ATM-0130, Ammonium Analysis in L-Histidine Monohydrochloride Monohydrate Via UPLC with UV Detection
- 6.2. BSI-PRL-0820, Analytical Method Validation Protocol: Ammonium Analysis in Amino Acids Via Ultra High-Performance Liquid Chromatography (UPLC) with UV Detection
- 6.3. BSI-RPT-1906, Analytical Method Validation Report: Ammonium Analysis in L-Histidine Monohydrochloride Monohydrate Via Ultra High-Performance Liquid Chromatography (UPLC) with UV Detection
- 6.4. BSI-SOP-0019, Result Reporting
- 6.5. BSI-SOP-0069, Preparation of Samples for Outside Testing
- 6.6. BSI-SOP-0098, Balance SOP
- 6.7. BSI-SOP-0126, Laboratory Notebooks
- 6.8. BSI-SOP-0134, Pipette SOP
- 6.9. BSI-SOP-0140, Standardization of Titrants
- 6.10. BSI-SOP-0143, Metrohm Titrando 907 Auto-Titrator SOP
- 6.11. BSI-SOP-0242, Bangor Portable Turbidimeter Operation and Calibration
- 6.12. BSI-SOP-0255, XL200 pH/mV/Conductivity Meter SOP
- 6.13. BSI-SOP-0490, MCP 5300 Polarimeter SOP
- 6.14. BSI-SOP-0668, OPI-180 OD Handheld Colorimeter SOP
- 6.15. Current USP/ NF
- 6.16. *EP Compendia*
- 6.17. ICH Q1A

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

6.18. *JP Compendia***7. ANALYTICAL PROCEDURES:****7.1. AMMONIUM (EP/JP) :**

7.1.1. Primary Analysis:

7.1.1.1. Refer to BSI-ATM-0130, Ammonium Analysis in L-Histidine Monohydrochloride Monohydrate via UPLC with UV Detection for sample preparation and analysis.

7.1.2. Alternative Analysis:

7.1.2.1. Refer to Section 7.15, Ninhydrin-Positive Substances.

7.2. APPEARANCE AND COLOR :

7.2.1. Place ~10 grams of sample into a clean, dry, glass beaker.

7.2.2. In an area with sufficient lighting, viewed the sample from all sides.

7.2.3. The sample should be white or colorless in color and characteristic of crystals or a crystalline powder.

7.2.4. If the appearance and color result is unable to be definitively determined visually, the sample may be analyzed using the Colorimeter. Refer to BSI-SOP-0668, OPI-OD Handheld Colorimeter SOP

7.2.5. Any non-conformance will be reported to the Laboratory Manager or designee, immediately.

7.3. APPEARANCE OF SOLUTION (EP) :7.3.1. Test Solution:

7.3.1.1. Dissolve 2.5 grams of sample in Purified Water, dilute to 50mL with Purified Water, and mix thoroughly.

7.3.2. Clarity:

7.3.2.1. Analyze the *Test Solution* for turbidity using a calibrated turbidimeter.

7.3.2.2. Acceptance Criteria:

7.3.2.2.1. The turbidity result may not exceed 3NTU to report as clear.

7.3.3. Color:

7.3.3.1. In an area with sufficient lighting, compare the color of the *Test Solution* to *EP Color Reference Solution BY₆*.

7.3.3.2. Acceptance Criteria:

7.3.3.2.1. The color of the *Test Solution* may not be more intense than the color in *EP Color Reference Solution BY₆* to report as colorless.

7.3.4. To report as Pass Test, the solution must be both Clear and Colorless.

7.4. ASSAY :

7.4.1. Perform a daily check or standardization of 0.1N Sodium Hydroxide as per Standardization of Titrants.

7.4.2. Accurately weigh 0.5 grams of sample into a suitable beaker.

7.4.3. Dissolve in 50mL of Purified Water.

7.4.4. Titrate to a potentiometric endpoint with 0.1N Sodium Hydroxide.

7.4.5. Calculate the %L-Histidine Monohydrochloride (as-is) assay value obtained from the Metrohm Tiamo software:

$$\% \text{L-Histidine Monohydrochloride (as-is)} = \frac{(\text{Sample Titrant Volume (mL)})(\text{Titrant Normality})(19.16)}{\text{Sample Weight (g)}}$$

7.4.6. Calculate the % L-Histidine Monohydrochloride Monohydrate (Dried Basis) and %L-Histidine Monohydrochloride Monohydrate (Anhydrous Basis) as needed using the following equations:

7.4.6.1. %L-Histidine Monohydrochloride Monohydrate (Dried Basis) calculation:

$$\%L\text{-Histidine Monohydrochloride Monohydrate (Dried Basis)} = \frac{\%L\text{-Histidine Monohydrochloride (as-is)} \times 100}{100 - LOD (\%)}$$

7.4.6.2. %L-Histidine Monohydrochloride Monohydrate (Anhydrous Basis) calculation:

$$\%L\text{-Histidine Monohydrochloride Monohydrate (Anhydrous Basis)} = \frac{\%L\text{-Histidine Monohydrochloride (as-is)} \times 100}{100 - Karl\ Fischer\ Value (\%)}$$

7.5. **BIOBURDEN** :

7.5.1. Package and send NLT 35 grams to Mary Paul Laboratories for TAMC/TYMC.

7.5.2. MPL Suitability Number to reference on ARF: 23P4414a

7.5.2.1. TAMC suitability <100CFU/g

7.5.2.2. TYMC suitability <10CFU/g

7.5.2.3. Bioburden Result Reported from TAMC result.

7.6. **CLARITY AND COLOR OF SOLUTION (JP)** :

7.6.1. Dissolve 1.0g of sample in 10mL of purified water. The solution must be clear and colorless.

7.7. **ENDOTOXIN** :

7.7.1. Accurately weigh 0.10 grams of sample into a sterile tube.

7.7.2. Dissolve in 5mL of LAL reagent water.

7.7.3. Add 0.4mL of 1M NaOH.

7.7.4. Dilute to 10mL with LAL reagent water for a final concentration of 0.010 g/mL.

7.7.5. Refer to Endosafe NexGen-PTS Endotoxin Reader SOP for instrument analysis.

7.8. **IDENTIFICATION (IR): C (EP), 1 (JP)** :

7.8.1. Follow the Spectrum Two UATR SOP for sample preparation and analysis.

7.8.2. Analyze sample as-is.

7.9. **IDENTIFICATION 2, CHLORIDE (JP)** :

7.9.1. Prepare a 1 in 10 sample solution by weighing 1 gram of sample and dissolving in 10mL of purified water. Mix thoroughly.

7.9.2. To the sample solution, add 0.2mL of 0.1N Silver Nitrate.

7.9.3. A white, curdy precipitate that is insoluble after the addition of 1mL of dilute Nitric Acid is produced.

7.9.4. If no precipitate is produced, notify Laboratory Management, or designee.

7.9.5. Add 4 mL of Ammonia TS. The precipitate should dissolve after minor agitation.

7.10. **IDENTIFICATION E (EP)** :

7.10.1. Weigh 0.1g of sample, transfer to a suitable beaker and dissolve in 7mL of purified water.

7.10.2. Add 3mL of 200g/L sodium hydroxide.

7.10.3. In a separate beaker, weigh 50mg of Sulfanilic acid and dissolve in a mixture of 0.1mL Hydrochloric Acid and 10mL of purified water. Add 0.1mL of sodium nitrite solution.

7.10.4. Add the Sulfanilic acid solution to the sample solution and mix.

7.10.5. An orange-red color must develop in order to report as passes test

7.11. **IDENTIFICATION F (EP)** :

7.11.1. Weigh 20mg of sample into a suitable vessel and dissolve in 2mL of purified water.

7.11.2. Acidify with dilute nitric acid.

7.11.3. Add 0.4mL of 0.1N Silver Nitrate.

7.11.3.1. A curdled, white precipitate is formed.

7.11.4. Centrifuge and wash the precipitate with three, 1mL quantities of purified water.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.11.4.1. Carry out rapidly in subdued light, disregarding the fact that the supernatant may not become clear.
- 7.11.5. To the precipitate, add 2mL of purified water and 1.5mL of ammonia R.
- 7.11.5.1. The precipitate dissolves easily with the possible exception of a few large particles which dissolve slowly in order to report as passes test.

7.12. HEAVY METALS (JP) :

- 7.12.1. Refer to BSI-ATM-0076, Analytical Method Determination of ICH Q3D Elemental Impurities (Class 1, 2A, 2B, 3 & 4) By Inductively Coupled Plasma Mass Spectrometry (ICP-MS) in L-Histidine Monohydrochloride Monohydrate for sample preparation and analysis.

7.13. IRON (EP/JP) :

- 7.13.1. Refer to BSI-ATM-0076, Analytical Method Determination of ICH Q3D Elemental Impurities (Class 1, 2A, 2B, 3 & 4) By Inductively Coupled Plasma Mass Spectrometry (ICP-MS) in L-Histidine Monohydrochloride Monohydrate for sample preparation and analysis.

7.14. LOSS ON DRYING (EP) :

- 7.14.1. Dry an LOD vial in an oven at 150°C for 30 minutes.
- 7.14.2. Cool for 15 minutes in a desiccator, weigh the LOD vial and record the weight.
 - 7.14.2.1. Note: If the substance to be tested is in the form of large crystals, reduce the particle size to about 2 mm by quickly crushing before weighing.
- 7.14.3. Transfer 1.000g of the sample to the LOD vial and accurately weigh the vial and contents. By gentle sideways shaking, distribute the sample as evenly as possible in the LOD vial to a depth of about 5mm.
- 7.14.4. Place the LOD vial containing the sample into the oven to dry at 150°C to a constant weight (NMT 0.5mg between weighings). Note: It is recommended to dry overnight to minimize the amount of times to reweigh the sample to ensure constant weight.
- 7.14.5. Remove the LOD vial from the oven and allow it to cool in the desiccator for 15 minutes.
- 7.14.6. Reweigh the LOD vial and sample and return to the oven for an additional 1 hour minimum.
- 7.14.7. Remove the LOD vial and sample from the oven and allow it to cool in the desiccator for 15 minutes.
- 7.14.8. Reweigh the LOD vial and sample and verify the weight has changed NMT 0.5mg between the last weighing.
- 7.14.9. Repeat steps 7.14.6-7.14.7 until a constant weight is achieved.
- 7.14.10. Retain the dried sample for dried analyses, if necessary.
- 7.14.11. Calculate the % LOD as follows:

$$\%LOD = \frac{\text{Initial Sample Weight (g)} - \text{Final Sample Weight (g)}}{\text{Initial Sample Weight (g)}} \times 100$$

7.15. NINHYDRIN-POSITIVE SUBSTANCES (EP) :

- 7.15.1. Package and send NLT 10g to New Jersey Labs for Ninhydrin-Positive substances.
 - 7.15.1.1. Method Reference on ARF: VER-21-045

7.16. OPTICAL ROTATION, SPECIFIC ROTATION @ 20°C; IDENTIFICATION A (EP) :

- 7.16.1. Sample Preparation:
 - 7.16.1.1. Weigh and transfer 11 grams of sample to a 100mL volumetric flask, dissolve in 48mL of Hydrochloric Acid R1, fill to volume with purified water, and mix well.
- 7.16.2. Optical Zero Reference:
 - 7.16.2.1. Pipette 48mL of Hydrochloric Acid R1 into a 100mL volumetric flask, fill to volume with purified water, and mix well.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.16.3. The following information will be required to be entered into the software:
 7.16.3.1. Volume (dryness) (mL), Mass (dryness) (g), Drying Loss (%)
 7.16.3.1.1. Utilize the sample's Loss on Drying result for Drying Loss (%).
- 7.16.4. Analysis: Perform at 20°C.
- 7.16.5. Refer to the MCP 5300 Polarimeter SOP for instrument analysis.

7.17. OPTICAL ROTATION, SPECIFIC ROTATION @ 20°C (JP) :

- 7.17.1. Sample Preparation:
 7.17.1.1. Weigh and transfer 11 grams of sample to a 100mL volumetric flask, dissolve in 6N Hydrochloric Acid, fill to volume with 6N Hydrochloric Acid, and mix well.
- 7.17.2. Optical Zero Reference:
 7.17.2.1. 6N Hydrochloric Acid
- 7.17.3. The following information will be required to be entered into the software:
 7.17.3.1. Volume (dryness) (mL), Mass (dryness) (g), Drying Loss (%)
 7.17.3.1.1. Utilize the sample's Karl Fischer Water Content result for Drying Loss (%).
- 7.17.4. Analysis: Perform at 20°C
- 7.17.5. Refer to the MCP 5300 Polarimeter SOP for instrument analysis.

7.18. pH (5% SOLUTION); IDENTIFICATION B (EP) :

- 7.18.1. Weigh 2.5g of sample and transfer to a beaker,
 7.18.2. Dissolve in 50mL of purified water and mix well.
 7.18.3. Follow the appropriate SOP for pH calibration and measurement.

7.19. pH (10% SOLUTION) (JP) :

- 7.19.1. Weigh 10.0g of sample and transfer to a suitable beaker.
 7.19.2. Dissolve in 100mL of purified water and mix well.
 7.19.3. Follow the appropriate SOP for pH calibration and measurement.

7.20. RELATED SUBSTANCES; IDENTIFICATION D (EP) :

- 7.20.1. Solution Preparation:
 7.20.1.1. Note: All solutions may be scaled as needed.
 7.20.1.2. Developing Solvent System (20:20:60 Glacial Acetic Acid: Purified Water: Butanol): Pipette 20mL of purified water into a 100mL volumetric flask. Add 20mL of glacial acetic acid. QS to the final volume with butanol.
 7.20.1.3. Spray Reagent: Ninhydrin Solution
 7.20.1.4. Sample Solution: Accurately weigh 10mg of sample and transfer to a 50mL volumetric flask. Dissolve and dilute to volume with purified water. Mix well.
 7.20.1.5. Standard Solution: Accurately weigh 10mg of L-Histidine Monohydrochloride Monohydrate CRS and transfer to a 50mL volumetric flask. Dissolve and dilute to volume with purified water.
- 7.20.2. Chromatographic System:

Table 1: Chromatographic System

Parameter	Setting
Mode	Thin-Layer Chromatography (TLC)
Adsorbent	0.25mm Layer of Chromatographic Silica Gel Mixture
Application Volume	5µL
Developing Solvent System	20:20:60 V/V/V Glacial Acetic Acid: Purified Water: Butanol
Spray Reagent	Ninhydrin R

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

7.20.3. Analysis:

- 7.20.3.1. Spot 5µL of *Sample Solution* and *Standard Solution* onto the TLC plate.
- 7.20.3.2. Place the plate in the chamber, ensuring the spots or bands are above the surface of the mobile phase.
- 7.20.3.3. Allow the mobile phase to ascent the plate until the solvent front has traveled two thirds of the length of the plate.
- 7.20.3.4. Remove the plate, mark the solvent front with a pencil.
- 7.20.3.5. Dry the plate at room temperature with air.
- 7.20.3.6. Spray with *Spray Reagent* and heat the plate at 105°C for 15 minutes.
- 7.20.3.7. Examine the plate under white light.
- 7.20.3.8. Acceptance Criteria: The principal spot for the test solution is similar in position, color, and size to the principal spot for the standard solution.

7.21. **RELATED SUBSTANCES (JP)** :

7.21.1. Solution Preparation:

- 7.21.1.1. Note: All solutions may be scaled as needed.
- 7.21.1.2. Developing Solvent System (67:33 1-Propanol: Ammonia Water): Add 67mL of 1-Propanol to 33mL of Ammonia Water.
- 7.21.1.3. Spray Reagent: Weigh 1 gram of ninhydrin and transfer to a 100mL volumetric flask. Dissolve and dilute to volume with a 97:3 mixture of methanol: acetic acid.
- 7.21.1.4. Sample Solution: Accurately weigh 0.50g of sample and transfer to a 50mL volumetric flask. Dissolve and dilute to volume with purified water. Mix well.
- 7.21.1.5. Standard Solution: Accurately weigh 0.10g of L-Histidine Monohydrochloride Monohydrate CRS and transfer to a 10mL volumetric flask.
 - 7.21.1.5.1. Pipette 1mL of the resulting solution into a 10mL volumetric flask and dilute to volume with purified water.
 - 7.21.1.5.2. Pipette 1mL of the second dilution solution into a 50mL volumetric flask and dilute to volume with purified water.

7.21.2. Chromatographic System:

Table 2: Chromatographic System

Parameter	Setting
Mode	Thin-Layer Chromatography (TLC)
Adsorbent	0.25mm Layer of Chromatographic Silica Gel Mixture
Application Volume	5µL
Developing Solvent System	67:33 1-Propanol: Ammonia Water
Spray Reagent	Ninhydrin in 97:3 Methanol: Acetic Acid

7.21.3. Analysis:

- 7.21.3.1. Spot 5µL of *Sample Solution* and *Standard Solution* onto the TLC plate.
- 7.21.3.2. Place the plate in the chamber, ensuring the spots or bands are above the surface of the mobile phase.
- 7.21.3.3. Allow the mobile phase to ascend the plate until the solvent front has traveled approximately 10cm.
- 7.21.3.4. Remove the plate, mark the solvent front with a pencil.
- 7.21.3.5. Dry the plate at 80°C for 30 minutes.
- 7.21.3.6. Spray with *Spray Reagent* and heat the plate at 80°C for 10 minutes.
- 7.21.3.7. Examine the plate under white light.
- 7.21.3.8. Acceptance Criteria: Any secondary spot of the sample solution is not more intense than the primary spot of the standard solution.

7.22. SULFATED ASH (EP), RESIDUE ON IGNITION (JP) :

- 7.22.1. Turn on the muffle furnace and allow it to stabilize at 600°C. Follow the muffle furnace calibration procedure for operation of the furnace.
- 7.22.2. Inspect a quartz crucible for cracks, chips, and discoloration.
- 7.22.3. Utilize forceps to insert and remove the crucible from the furnace.
- 7.22.4. Ignite a quartz crucible at $600 \pm 50^\circ\text{C}$ for 30 minutes. Cool in a desiccator for 1.5 hours and weigh using an analytical balance.
- 7.22.5. Weigh 1.0 grams of sample in the previously ignited quartz crucible. Moisten the sample with 1.0mL of concentrated Sulfuric Acid.
- 7.22.6. Volatilize the sample until the sample is thoroughly charred and white fumes are no longer evolved. Heat the sample slowly, so that the sample does not boil over and sample is not lost.
 - 7.22.6.1. The rate of heating should be such that from $\frac{1}{2}$ to 1 hour is required to volatilize the sample.
 - 7.22.6.2. Continue to heat the sample until all excess sulfuric acid has been volatilized.
- 7.22.7. Ignite in the muffle furnace at $600 \pm 50^\circ\text{C}$ for 15 minutes or until all carbon has been removed.
- 7.22.8. Cool in a desiccator for the same amount of time employed in the preparation of the crucible and weigh on an analytical balance.
- 7.22.9. Calculate the %ROI as follows:

$$\%ROI = \frac{\text{Residue Weight (g)}}{\text{Sample Weight (g)}} \times 100$$

- 7.22.10. If the amount of residue exceeds the limit specified, repeat the moistening with sulfuric acid using up to 1mL, heat to char, then ignite at $600 \pm 50^\circ\text{C}$ for 30 minutes until two consecutive weighings of the residue do not differ by more than 0.0005g or until the specification is met.

7.23. SULFATES (EP) :

- 7.23.1. Solution S: Weigh 2.5g of sample and transfer to a 50mL volumetric flask. Dissolve and dilute to volume with purified water.
- 7.23.2. Test Solution: Dilute 10mL of Solution S to 15mL with purified water in a suitable beaker.
- 7.23.3. Sulfate Standard Solution (10ppm SO_4) R1: Immediately before use, weigh 0.181g of dipotassium sulfate into a 100mL volumetric flask and dilute to volume with reagent alcohol (30% V/V). Pipette 1mL of the resulting solution into a 100mL volumetric flask and dilute to volume with reagent alcohol (30% V/V).
- 7.23.4. Sulfate Standard Solution (10ppm SO_4): Immediately before use, weigh 0.181g of dipotassium sulfate into a 100mL volumetric flask and dilute to volume with purified water. Pipette 1mL of the resulting solution into a 100mL volumetric flask and dilute to volume with purified water.
- 7.23.5. Procedure: In a beaker, add 3mL of 250 g/L Barium Chloride to 4.5mL Sulfate Standard Solution (10ppm SO_4) R1. Shake and allow to stand for 1 minute. Pipette 2.5mL of the resulting solution into the 15mL of test solution and add 0.5mL of acetic acid. Prepare a standard in the same manner using 15mL of Sulfate Standard Solution (10ppm SO_4) previously prepared instead of the test solution. Allow solutions to sit for 5 minutes. Any opalescence in the test solution is not more intense than that of the standard in order to report as <300 ppm.

7.24. **SULFATES (JP)** _____ :

7.24.1. Standard Preparation:

7.24.1.1. Pipette 0.35mL of 0.005M Sulfuric Acid into a 50mL Nessler Color Comparison Tube. Add ~40mL of purified water.

7.24.1.2. Add 1 mL of dilute hydrochloric acid and dilute to 50mL with purified water.

7.24.2. Sample Preparation:

7.24.2.1. Weigh 0.6 grams of sample and transfer to a 50mL Nessler Color Comparison Tube. Dissolve in ~40mL of purified water.

7.24.2.2. Add 1mL of dilute hydrochloric acid and dilute to 50mL with purified water.

7.24.3. Procedure:

7.24.3.1. To both the standard and the sample, add 2mL of barium chloride TS.

7.24.3.2. Mix well and allow to stand for 10 minutes utilizing a calibrated timer.

7.24.3.3. Compare the turbidity produced in both solutions against a black background by viewing downward or transversely.

7.24.3.4. Any turbidity produced in the sample solution should not exceed that produced by the standard in order to report as <0.028%.

7.25. **WATER** _____ :

7.25.1. Standardize Composite 5 as per Standardization of Titrants.

7.25.2. Grind the sample in a dry mortar into a fine powder utilizing a pestle.

7.25.3. Immediately weigh 0.3 g of sample into the glass weighing spoon and tare it.

7.25.4. Transfer the sample to the KF vessel by removing the rubber septum and adding the sample into the titration vessel.

7.25.4.1. Do not leave the rubber septum open for long periods of time as this will allow moisture to enter the titration vessel.

7.25.5. Return the weighing spoon to the balance, making sure not to lose any sample that was left behind. Once the weight stabilizes, record the weight into the Tiamo software.

7.25.6. Check to make sure there is no residual sample stuck to the sides of the titration vessel.

7.25.6.1. If there is any sample stuck to the side, stop the stir bead from spinning before swirling the vessel to rinse the sides.

7.25.7. Once the method begins, check to ensure the sample is fully dissolved before the titration begins (i.e., before the stir command completes).

7.25.8. The moisture content will then be determined by the Metrohm Titrando 907.

$$\%Moisture = \frac{(mL \text{ of Composite 5}) \left(\frac{mg}{mL} \text{ of Composite 5} \right) (0.1)}{Sample \text{ weight (g)}}$$

8. COMPENDIAL DIFFERENTIATIONS:

8.1. Compendial Analyses Not Harmonized

Table 3: Compendial Analyses Not Harmonized

EP Compendia	JP Compendia
Appearance of Solution	Clarity and Color of Solution
Identification, Specific Optical Rotation (EP-A)	Identification 2, Chloride
Identification B, pH	Optical Rotation
Identification D	pH
Identification E	Related Substances
Identification F	Sulfates
Loss on Drying	
Sulfates	

8.2. One Method Utilized

Table 4: One Method Name

Analysis Name
Identification, IR
Residue on Ignition/Sulfated Ash

8.3. In-House Validated Methods in accordance with USP General Chapter requirements:

Table 5: In-House Validated Methods

Analysis Name
Ammonium and Ninhydrin-positive substances
Assay (different calculations to report dried substance and anhydrous basis)
Endotoxin
Trace Metals to report: Heavy Metals and Iron
Water

8.4. In-House Methods for Product Quality Description

Table 6: In-House Methods

Analysis Name
Appearance and Color

8.5. Outside Approved Laboratory Testing (if required)

Table 7: Outside Approved Laboratory Testing

Analysis Name	
Bioburden (TAMC/TYMC) Sample size required to be sent to MPL:~35grams	MPL Suitability #23P4414a
Ammonium and Ninhydrin-positive substances Sample size required to be sent to NJ labs: ~10grams (Outside lab can be utilized as alternate method if required)	Method Reference: VER-21-045

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.