

UREA 6M SOLUTION

GMP, EXCIPIENT

CAS #: 57-13-6

Formula: $\text{CH}_4\text{N}_2\text{O}$

M.W.: 60.06 g/mol

UREA-3120**Bio EXCIPIENT GRADE**

ANALYSIS	SPECIFICATIONS
Appearance and Color	Colorless Solution
Identification, IR	Conforms to Reference Standard
Molarity	5.8 – 6.2 M
pH, 25°C	7 – 10
Trace Metals	Arsenic (As) < = 5 ppm Copper (Cu) < = 5 ppm Iron (Fe) < = 5 ppm Lead (Pb) < = 5 ppm

General Product Overview

Urea is used as a protein denaturant with low UV absorptivity. In addition to increasing the solubility of hydrophobic molecules, Urea denatures proteins by unfolding them and altering their three-dimensional structure. BioSpectra manufactures GMP, multi-compendial Urea crystal which is combined with water for injection to manufacture the finished 6M Urea Solution.

Industry Application

Suitable for use as a cGMP chemical in pharmaceutical manufacturing processes and products.

Key Product Features

- The manufacturing of Urea 6M Solution is performed at BioSpectra's Bangor, PA facility.
- Appears as a colorless solution
- Manufactured in accordance with ICH Q7
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: Carbamide Solution; Carbonyl Diamide Solution

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Unless otherwise noted on the Shelf-Life Statement and CoA, this product has a 2-year retest date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

200L drums and 1,000L – 1,200L totes

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.