

## TRIS HYDROCHLORIDE

## PROCESS FLOW DIAGRAM - BANGOR

Mother Liquor Creation or Verification/Adjustment

Tris and Tris HCl Charged into Mother Liquor

Hydrochloric Gas Addition to Correct pH (Exothermic Reaction)



- 200 - 1000000

Recrystallization

Wet Crystal Separation via Centrifuge

Drying to Desired Moisture Through Fluid Bed Dryer

Final Packaging

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.