

RESIDUAL SOLVENT STATEMENT

L-Glutamine

L-Glutamine, Bio Qualified Grade is manufactured using the solvent Ethyl Acetate, and is stated to comply with the ICH Q3C Guideline for Residual Solvents required concentration limit of 5000ppm for a Class 3 residual solvent. BioSpectra has evaluated the material supply through the Supplier Qualification Program. BioSpectra has additionally analyzed L-Glutamine, Bio Qualified Grade for residual solvents with results meeting the specifications.

Current Product Number
LGLM-6201

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist