

# HEPES, FREE ACID

## LBLE, GMP

CAS #: 7365-45-9

Formula:  $C_8H_{18}N_2O_4S$

F.W.: 238.30 g/mol

### HEPE-4250

### BIO PHARMA GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (0.1M)	250 nm	< = 0.0500 a.u.
	260 nm	< = 0.0500 a.u.
	280 nm	< = 0.0800 a.u.
Absorbance (0.05M)	250 nm	< = 0.0500 a.u.
	260 nm	< = 0.0500 a.u.
	280 nm	< = 0.0800 a.u.
Appearance and Color		White Crystals
Assay, Dried Basis		> = 99.5%
Chloride		< = 50 ppm
Endotoxin		< = 5 EU/g
Enzymes	DNase	None Detected
	RNase	None Detected
	Protease	None Detected
Identification, IR		Conforms to Reference Standard
Insoluble Matter		< = 0.01%
Loss on Drying		< = 0.5%
Microbial Content	TAMC	< = 100 CFU/g
	TYMC	< = 100 CFU/g
pH (5%)		5.0 – 6.5
pK <sub>a</sub>		7.45 – 7.65
Residue on Ignition		< = 0.1%
Solubility	Solubility (5%)	Passes Test
	Solubility (0.05M)	Passes Test
Sulfate		< = 50 ppm
Trace Metals	Arsenic (As)	< = 5 ppm
	Aluminum (Al)	< = 5 ppm
	Bismuth (Bi)	< = 5 ppm
	Calcium (Ca)	< = 10 ppm
	Copper (Cu)	< = 5 ppm
	Iron (Fe)	< = 5 ppm
	Lead (Pb)	< = 1 ppm
	Lithium (Li)	< = 5 ppm
	Molybdenum (Mo)	< = 5 ppm
	Nickel (Ni)	< = 5 ppm
Water, KF	Potassium (K)	< = 50 ppm
		< = 0.1%

### General Product Overview

This is a high purity, reagent grade product, purified under cGMP conditions. HEPES, Free Acid is a zwitterionic buffer used to maintain pH of media used in cell cultures. It is one of Good's buffers that has a pK<sub>a</sub> value similar to its pH value, making it an ideal buffer for pH maintenance. HEPES is a Good's buffer because it has low UV absorptivity, minimal reactivity, stable pH and is soluble in water.

### Industry Application

Suitable for use in biological and biotech chemical process applications from R&D through scale cGMP production.

### Key Product Features

- Appears as a white crystalline product
- Manufactured in accordance with IPEC
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: N-(2-Hydroxyethyl) Piperazine-N'-2-Ethanesulfonic Acid; 4-(2-Hydroxyethyl) Piperazine-1-Ethanesulfonic Acid

### Storage and Shipping Conditions

Refer to SDS.

### Standard Shelf-Life Policy

Unless otherwise noted on the Shelf-Life Statement and CoA, this product has a 2-year retest date supported by a 3-year ICH Q1 Stability Study (if one is completed).

### Package Sizes

100g, 500g, 1kg, 5kg, 10kg, 25kg, 50kg

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*